

Safety Data Sheet

According to EC Directive 91/155/EEC

Date of issue: 30.04.2004
Supersedes edition of 28.01.2000

1. Identification of the substance/preparation and of the company/undertaking

Identification of the product

Catalogue No.: 114848
Product name: Phosphate Test Method: photometric, PMB
also refill pack for art. 1.14846 0.010 - 5.00 mg/l PO₄-P 420
Tests
0.03 - 15.3 mg/l PO₄³⁻
0.02 - 11.46 mg/l P₂O₅
Spectroquant®
P-1A

Use of the substance/preparation

Reagent for analysis

Company/undertaking identification

Company: Merck KGaA * 64271 Darmstadt * Germany * Tel: +49 (0)6151/72-0
Emergency telephone No.: +49 (0)6151/72112 * Fax: +49 (0)6151/72-7780

2. Composition/information on ingredients

Preparation contains inorganic compounds.

Hazardous ingredients:

Name according to EC Directives:

CAS-No.	EC No.	EC-Index-No.	Classification	Content:
Sulphuric acid				
7664-93-9	231-639-5	016-020-00-8	C; R35	≥ 25 - < 50 %

(Full text of R-Phrases in heading 16)

3. Hazards identification

Causes severe burns.

4. First aid measures

After inhalation: fresh air. Call in physician.
After skin contact: wash off with plenty of water. Dab with polyethylene glycol 400.
Immediately remove contaminated clothing.
After eye contact: rinse out with plenty of water for at least 10 minutes with the eyelid held wide open. Immediately call in ophthalmologist.
After swallowing: make victim drink plenty of water (if necessary several litres), avoid vomiting (risk of perforation!). Immediately call in physician. Do not attempt to neutralize.

Merck Safety Data Sheet

According to EC Directive 91/155/EEC

Catalogue No.: 114848
Product name: Phosphate Test Method: photometric, PMB
also refill pack for art. 1.14846 0.010 - 5.00 mg/l PO₄-P 420 Tests
0.03 - 15.3 mg/l PO₄³⁻
0.02 - 11.46 mg/l P₂O₅
Spectroquant®
P-1A

5. Fire-fighting measures

Suitable extinguishing media:

In adaption to materials stored in the immediate neighbourhood.

Special risks:

Non-combustible. Ambient fire may liberate hazardous vapours. The following may develop in event of fire: sulfur oxides.

Special protective equipment for fire fighting:

Do not stay in dangerous zone without self-contained breathing apparatus. In order to avoid contact with skin, keep a safety distance and wear suitable protective clothing.

Other information:

Contain escaping vapours with water. Prevent fire-fighting water from entering surface water or groundwater.

6. Accidental release measures

Person-related precautionary measures:

Do not inhale vapours/aerosols. Avoid substance contact. Ensure supply of fresh air in enclosed rooms.

Environmental-protection measures:

Do not allow to enter sewerage system.

Procedures for cleaning / absorption:

Take up with liquid-absorbent and neutralizing material (e.g. Chemizorb® H⁺, Art. No. 101595). Forward for disposal. Clean up affected area.

7. Handling and storage

Handling:

No further requirements.

Storage:

Tightly closed. At +15°C to +25°C.

The data apply to the entire pack.

8. Exposure controls/personal protection

Specific control parameter

Merck Safety Data Sheet

According to EC Directive 91/155/EEC

Catalogue No.: 114848
Product name: Phosphate Test Method: photometric, PMB
also refill pack for art. 1.14846 0.010 - 5.00 mg/l PO₄-P 420 Tests
0.03 - 15.3 mg/l PO₄³⁻
0.02 - 11.46 mg/l P₂O₅
Spectroquant®
P-1A

TRGS 900

Name	Sulfuric acid
Kind of use	others
Value	0.1 mg/m ³ Inhalable fraction. The following regulation applies to the peak limit: short-term exposure duration: max. 15 minutes as mean value, frequency per shift: 4, time interval: min. 1 hour
Peak limit	1 Concentration must not exceed limit concentration.
Embryotoxic	Y Substances with which no foetotoxic risk is to be expected when observing the maximum allowable concentration (MAC Germany) and the biological tolerance value at the workplace (BAT Germany).

Personal protective equipment:

Protective clothing should be selected specifically for the working place, depending on concentration and quantity of the hazardous substances handled. The resistance of the protective clothing to chemicals should be ascertained with the respective supplier.

Respiratory protection: required when vapours/aerosols are generated.

Eye protection: required

Hand protection:

In full contact:

Glove material: natural latex
Layer thickness: 0.6 mm
Breakthrough time: > 480 Min.

In splash contact:

Glove material: natural latex
Layer thickness: 0.6 mm
Breakthrough time: > 480 Min.

The protective gloves to be used must comply with the specifications of EC directive 89/686/EEC and the resultant standard EN374, for example KCL 706 Lapren® (full contact), 706 Lapren® (splash contact). The breakthrough times stated above were determined by KCL in laboratory tests acc. to EN374 with samples of the recommended glove types.

This recommendation applies only to the product stated in the safety data sheet and supplied by us as well as to the purpose specified by us. When dissolving in or mixing with other substances and under conditions deviating from those stated in EN374 please contact the supplier of CE-approved gloves (e.g. KCL GmbH, D-36124 Eichenzell, Internet: www.kcl.de).

Industrial hygiene:

Change contaminated clothing and immerse in water. Apply skin-protective barrier cream. Wash hands and face after working with substance.

Merck Safety Data Sheet

According to EC Directive 91/155/EEC

Catalogue No.: 114848
Product name: Phosphate Test Method: photometric, PMB
also refill pack for art. 1.14846 0.010 - 5.00 mg/l PO₄-P 420 Tests
0.03 - 15.3 mg/l PO₄³⁻
0.02 - 11.46 mg/l P₂O₅
Spectroquant®
P-1A

9. Physical and chemical properties

Form:	liquid	
Colour:	colourless	
Odour:	odourless	
pH value	(20 °C)	acid
Melting point		not available
Boiling point		not available
Ignition temperature		not available
Flash point		not available
Explosion limits	lower	not available
	upper	not available
Density	(20 °C)	~ 1.21 g/cm ³
Solubility in water	(20 °C)	soluble

10. Stability and reactivity

Conditions to be avoided

Heating.

Substances to be avoided

water, alkali metals, alkali compounds, ammonia, alkaline earth metals, alkaline earth compounds, alkalis, acids, metals, metal alloys, phosphorus, phosphorus oxides, hydrides, halogen-halogen compounds, oxyhalogenic compounds, permanganates, nitrates, carbides, combustible substances, organic solvents, acetylidene, nitriles, organic nitro compounds, anilines, peroxides, picrates, nitrides, lithium silicide.

Hazardous decomposition products

in the event of fire: See chapter 5.

Further information

hygroscopic; has a corrosive effect
incompatible with various metals, animal/vegetable tissues.

11. Toxicological information

Acute toxicity

Quantitative data on the toxicity of this product are not available.

Merck Safety Data Sheet

According to EC Directive 91/155/EEC

Catalogue No.: 114848
Product name: Phosphate Test Method: photometric, PMB
also refill pack for art. 1.14846 0.010 - 5.00 mg/l PO₄-P 420 Tests
0.03 - 15.3 mg/l PO₄³⁻
0.02 - 11.46 mg/l P₂O₅
Spectroquant®
P-1A

Further toxicological information

Property that must be anticipated on the basis from the components of the preparation:

After inhalation of aerosols: damage to the affected mucous membranes.

After skin contact: severe burns with formation of scabs.

After eye contact: burns, corneal lesions.

After swallowing: severe pain (risk of perforation!), nausea, vomiting and diarrhoea.

After a latency period of several weeks possibly pyloric stenosis.

Further data

Further hazardous properties cannot be excluded.

The product should be handled with the care usual when dealing with chemicals.

12. Ecological information

Ecotoxic effects:

Quantitative data on the ecological effect of this product are not available.

Further ecologic data:

The following applies to sulfuric acid: biological effects: harmful effect on aquatic organisms.

Harmful effect due to pH shift. Toxic effect on fish and algae. Caustic even in diluted form.

Does not cause biological oxygen deficit. Endangers drinking-water supplies if allowed to enter soil and/or waters in large quantities. Neutralization possible in waste water treatment plants.

Daphnia toxicity: Daphnia magna EC₅₀: 29 mg/l/24 h (calculated on the pure substance).

Do not allow to enter waters, waste water, or soil!

13. Disposal considerations

Product:

Chemicals must be disposed of in compliance with the respective national regulations. Under www.retrologistik.de you will find country- and substance-specific information as well as contact partners.

Packaging:

Merck product packaging must be disposed of in compliance with the country-specific regulations or must be passed to a packaging return system. Under www.retrologistik.de you will find special information on the respective national conditions as well as contact partners.

Merck Safety Data Sheet

According to EC Directive 91/155/EEC

Catalogue No.: 114848
Product name: Phosphate Test Method: photometric, PMB
also refill pack for art. 1.14846 0.010 - 5.00 mg/l PO₄-P 420 Tests
0.03 - 15.3 mg/l PO₄³⁻
0.02 - 11.46 mg/l P₂O₅
Spectroquant®
P-1A

14. Transport information

Road & Rail ADR, RID
UN 3316 CHEMIE-TESTSATZ, 9, II
Inland waterway ADN, ADNR not tested
Sea IMDG-Code
UN 3316 CHEMICAL KIT, 9, II
Ems F-A S-P
Air CAO, PAX
CHEMICAL KIT, 9, UN 3316, II

The transport regulations are cited according to international regulations and in the form applicable in Germany. Possible national deviations in other countries are not considered.
THESE TRANSPORT DATA APPLY TO THE ENTIRE PACK !

15. Regulatory information

Labelling according to EC Directives

Symbol:	C	Corrosive
R-phrases:	35	Causes severe burns.
S-phrases:	26-36/37/39-45	In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. Wear suitable protective clothing, gloves and eye/face protection. In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).

contains: Sulphuric acid

German regulations

Water pollution class 1 (slightly polluting substance) VwVwS Anh. 4

Storage class VCI 8 B

Data sheet of the Chemical Professional Association M004 Irritant substances/corrosive substances
M051 Dangerous chemical substances

The employment restrictions for young workers in accordance with section 22 of the Youth Employment Protection Law (JArbSchG) are to be observed.

16. Other information

Text of any R phrases referred to under heading 2:

35 Causes severe burns.

Merck Safety Data Sheet

According to EC Directive 91/155/EEC

Catalogue No.: 114848
Product name: Phosphate Test Method: photometric, PMB
also refill pack for art. 1.14846 0.010 - 5.00 mg/l PO₄-P 420 Tests
0.03 - 15.3 mg/l PO₄³⁻
0.02 - 11.46 mg/l P₂O₅
Spectroquant®
P-1A

Reduced labelling on the container due to small quantity.

Reason for alteration

Chapter 2: change in ingredients data.
Chapter 14: transport information.

General update.

Contact for information:

HSSE-C/CI * Tel: +49 (0)6151/722775 * Fax: +49 (0)6151/726433 * e-mail:prodsafe@merck.de

The information contained herein is based on the present state of our knowledge. It characterizes the product with regard to the appropriate safety precautions. It does not represent a guarantee of the properties of the product.

Safety Data Sheet

According to EC Directive 91/155/EEC

Date of issue: 30.04.2004
Supersedes edition of 28.01.2000

1. Identification of the substance/preparation and of the company/undertaking

Identification of the product

Catalogue No.: 114848
Product name: Phosphate Test Method: photometric, PMB
also refill pack for art. 1.14846 0.010 - 5.00 mg/l PO₄-P 420
Tests
0.03 - 15.3 mg/l PO₄³⁻
0.02 - 11.46 mg/l P₂O₅
Spectroquant®
P-2A

Use of the substance/preparation

Reagent for analysis

Company/undertaking identification

Company: Merck KGaA * 64271 Darmstadt * Germany * Tel: +49 (0)6151/72-0
Emergency telephone No.: +49 (0)6151/72112 * Fax: +49 (0)6151/72-7780

2. Composition/information on ingredients

CAS-No.: 50-81-7
M: 176.13 g/mol
Formula Hill: C₆H₈O₆
EC-No.: 200-066-2

3. Hazards identification

No hazardous product as specified in Directive 67/548/EEC.

4. First aid measures

After inhalation: fresh air.
After skin contact: wash off with plenty of water. Remove contaminated clothing.
After eye contact: rinse out with plenty of water with the eyelid held wide open. Call in ophthalmologist if necessary.
After swallowing (large amounts): consult doctor if feeling unwell.

5. Fire-fighting measures

Suitable extinguishing media:
foam, powder, water.

Special risks:
Combustible. Development of hazardous combustion gases or vapours possible in the event of fire.

Special protective equipment for fire fighting:
Do not stay in dangerous zone without self-contained breathing apparatus.

Other information:
Prevent fire-fighting water from entering surface water or groundwater.

Merck Safety Data Sheet

According to EC Directive 91/155/EEC

Catalogue No.: 114848
Product name: Phosphate Test Method: photometric, PMB
also refill pack for art. 1.14846 0.010 - 5.00 mg/l PO₄-P 420 Tests
0.03 - 15.3 mg/l PO₄³⁻
0.02 - 11.46 mg/l P₂O₅
Spectroquant®
P-2A

6. Accidental release measures

Person-related precautionary measures:
Avoid generation of dusts; do not inhale dusts.

Environmental-protection measures:
Do not allow to enter sewerage system.

Procedures for cleaning / absorption:
Take up dry. Forward for disposal. Clean up affected area.

7. Handling and storage

Handling:

Protect from light.

Storage:

Tightly closed. Dry. At +15°C to +25°C.

The data apply to the entire pack.

8. Exposure controls/personal protection

Personal protective equipment:

Protective clothing should be selected specifically for the working place, depending on concentration and quantity of the hazardous substances handled. The resistance of the protective clothing to chemicals should be ascertained with the respective supplier.

Respiratory protection: required when dusts are generated.

Eye protection: required

Merck Safety Data Sheet

According to EC Directive 91/155/EEC

Catalogue No.: 114848
Product name: Phosphate Test Method: photometric, PMB
also refill pack for art. 1.14846 0.010 - 5.00 mg/l PO₄-P 420 Tests
0.03 - 15.3 mg/l PO₄³⁻
0.02 - 11.46 mg/l P₂O₅
Spectroquant®
P-2A

Hand protection:

In full contact:
Glove material: nitrile rubber
Layer thickness: 0.11 mm
Breakthrough time: > 480 Min.

In splash contact:
Glove material: nitrile rubber
Layer thickness: 0.11 mm
Breakthrough time: > 480 Min.

The protective gloves to be used must comply with the specifications of EC directive 89/686/EEC and the resultant standard EN374, for example KCL 740 Dermatril® (full contact), 740 Dermatril® (splash contact). The breakthrough times stated above were determined by KCL in laboratory tests acc. to EN374 with samples of the recommended glove types.

This recommendation applies only to the product stated in the safety data sheet and supplied by us as well as to the purpose specified by us. When dissolving in or mixing with other substances and under conditions deviating from those stated in EN374 please contact the supplier of CE-approved gloves (e.g. KCL GmbH, D-36124 Eichenzell, Internet: www.kcl.de).

Industrial hygiene:
Change contaminated clothing. Wash hands after working with substance.

9. Physical and chemical properties

Form:	solid
Colour:	white
Odour:	odourless
pH value at 50 g/l H ₂ O	(20 °C) 2.2-2.5
Melting point	190-192 °C (decomposition)
Boiling point	not applicable
Ignition temperature	380 °C
Flash point	not applicable
Explosion limits	lower not applicable upper not applicable
Vapour pressure	not applicable
Density	1.65 g/cm ³
Bulk density	~ 500-900 kg/m ³
Solubility in water	(24 °C) 330 g/l
Thermal decomposition	> 192 °C
log Pow:	-2.15

Merck Safety Data Sheet

According to EC Directive 91/155/EEC

Catalogue No.: 114848
Product name: Phosphate Test Method: photometric, PMB
also refill pack for art. 1.14846 0.010 - 5.00 mg/l PO₄-P 420 Tests
0.03 - 15.3 mg/l PO₄³⁻
0.02 - 11.46 mg/l P₂O₅
Spectroquant®
P-2A

10. Stability and reactivity

Conditions to be avoided

Heating.

Substances to be avoided

aluminium, copper, copper alloys, zinc, metal ions, oxidizing agent.

Hazardous decomposition products

no information available

Further information

light-sensitive; sensitive to air; sensitive to moisture.
reducing agents.

11. Toxicological information

Acute toxicity

LD₅₀ (oral, rat): 11900 mg/kg.

Specific symptoms in animal studies:

Eye irritation test (rabbit): Slight irritations.

Skin irritation test (rabbit): No irritation.

Further toxicological information

After eye contact: Slight irritations.

Further data

Naturally occurring substance.

The product should be handled with the care usual when dealing with chemicals.

Merck Safety Data Sheet

According to EC Directive 91/155/EEC

Catalogue No.: 114848
Product name: Phosphate Test Method: photometric, PMB
also refill pack for art. 1.14846 0.010 - 5.00 mg/l PO₄-P 420 Tests
0.03 - 15.3 mg/l PO₄³⁻
0.02 - 11.46 mg/l P₂O₅
Spectroquant®
P-2A

12. Ecological information

Biologic degradation:
Readily biodegradable.
When introduced properly, no impairments in the function of adapted biological waste-water-treatment plants are to be expected.

Behavior in environmental compartments:
Distribution: log p(o/w): -2.15.
No bioaccumulation is to be expected (log P(o/w) <1).

Ecotoxic effects:
Biological effects:
Fish toxicity:
Onchorhynchus mykiss LC₅₀: 1020 mg/l /96 h acid;
L.idus LC₅₀: 33000 mg/l /48 h (own results).

Further ecologic data:
BOD 48 % from TOD /5 d (closed bottle test);
BOD 65 % from TOD /28 d (closed bottle test).

No ecological problems are to be expected when the product is handled and used with due care and attention.

13. Disposal considerations

Product:

Chemicals must be disposed of in compliance with the respective national regulations. Under www.retrologistik.de you will find country- and substance-specific information as well as contact partners.

Packaging:

Merck product packaging must be disposed of in compliance with the country-specific regulations or must be passed to a packaging return system. Under www.retrologistik.de you will find special information on the respective national conditions as well as contact partners.

14. Transport information

Road & Rail ADR, RID
UN 3316 CHEMIE-TESTSATZ, 9, II
Inland waterway ADN, ADNR not tested
Sea IMDG-Code
UN 3316 CHEMICAL KIT, 9, II
Ems F-A S-P
Air CAO, PAX
CHEMICAL KIT, 9, UN 3316, II

The transport regulations are cited according to international regulations and in the form applicable in Germany. Possible national deviations in other countries are not considered.
THESE TRANSPORT DATA APPLY TO THE ENTIRE PACK !

Merck Safety Data Sheet

According to EC Directive 91/155/EEC

Catalogue No.: 114848
Product name: Phosphate Test Method: photometric, PMB
also refill pack for art. 1.14846 0.010 - 5.00 mg/l PO₄-P 420 Tests
0.03 - 15.3 mg/l PO₄³⁻
0.02 - 11.46 mg/l P₂O₅
Spectroquant®
P-2A

15. Regulatory information

Labelling according to EC Directives

Symbol: ---

R-phrases: ---

S-phrases: ---

German regulations

Water pollution class 1 (slightly polluting substance) VwVwS Anh. 2 KennNr. 737

Storage class VCI 10-13

16. Other information

Reason for alteration

Chapter 14: transport information.

General update.

Contact for information:

HSSE-C/CI * Tel: +49 (0)6151/722775 * Fax: +49 (0)6151/726433 * e-mail:prodsafe@merck.de

The information contained herein is based on the present state of our knowledge. It characterizes the product with regard to the appropriate safety precautions. It does not represent a guarantee of the properties of the product.