# Basic

# CM • TEB

## Pocket- • Powder- • Carat balance



Pocket balance

STANDARD

EXT



### Standard



Adjusting program (CAL): For quick setting of the balance's accuracy. External adjusting weight required.



**Battery operation:** The battery type is specified for each device.



TAR Tare: put back displayed value to "0" during weighing process, e.g. to add or remove weighing material from a container.



Small, practical and protected in a hard case cover

> Hard case cover as protection against pressure and dust. With integrated

Model	range	g Read- out		
KERN	Max	<b>d</b> g		
CM 60-2N	60	0,01		
CM 150-1N	150	0,1		
CM 320-1N	320	0,1		
CM 1K1N	1000	1		
	000	0.1		
TEB 200-1	200	0,1		

Overall dimensions WxDxH 85x130x25 mm. Net weight approx. 180 g. Backlit LCD display (only KERN CM-N) digit height 12 mm.

Powder balance

with grain-division (gn)

CM 500-GN

Now also for jewellers in carat, for reloaders in grain

Carat balance . Powder balance

STANDARD

EXT



Can be switched over from ct to g at the touch of a key

Ready for use Batteries size AAA (2x 1.5 V) included. AUTO-OFF

function after 4 min without a change of load to preserve the batteries, can be switched off

Test weight included Option: With DKD-Test weight-Calibration certificate, see page 91, KERN 960-477

#### Ready for use

Batteries size AAA (2x 1.5 V) included. AUTO-OFF function after 4 min without a change of load to preserve the batteries, can be switched off

# **Draft shield**

against air movements standard (both models)

Hard case cover as protection against pressure

Weighing pan included (both models)

Can be switched over from gn to g at the touch of a key

Price reduction

Weighing Read-Model out KERN ct | g ct CM 50-C2N 50 | 10 0,01 | 0,002 gn | g gn | g CM 500-GN1 500 | 32 0,1 | 0,01

Overall dimensions WxDxH 85x130x25 mm. Net weight approx. 200 g 13 and dust

# KERN - Your DKD-Calibration Partner

issues DKD calibration certificates for balances and test weights in its own accredited laboratories. It is internationally valid. KERN is a manufacturer of test

weights in the international classes of accuracy E1, E2, F1, F2, M1, M2, M3, from 1 mg up to 2000 kg DKD = German

Calibration Service



# Pour No.

## Your KERN specialist dealer



ADVANCED APPLIED TECHNOLOGIES Contact Us:

Irl Ph: 01 4523432 UK Ph: 08452 30 40 30 Web: www.carlstuart.com Email: info@carlstuart.com