Basic



Standard




Adjusting program (CAL):
For quick setting of the balance's accuracy. External adjusting weight required.



Suspended weighing: load support with hook on the underside of the balance.

Battery operation:



0.00

TAR Tare: put back displayed value to "0" during weighing process, e.g. to add or remove weighing

Option (additional price)

material from a container.



V 230 Mains adapter: 230 V/50Hz in standard version for Germany. On request GB, AUS or USA version.

EMB School balance





Chemistry school experiments

Weighing plate B

2 x 1.5 V, size AA (A). AUTO-OFF function after 3 min without a change of load to preserve the batteries, can be switched off



Mains adapter not standard, can be reordered at any time, KERN 440-902

15

Model	Weighing range	Read- out	Weighing plate		
KERN	9	g	mm		
EMB 200-2 =	200	0,01	A	9	
EMB 600-2 =	600	0,01	A		
EMB 220-1 A	220	0,1	В	9	
EMB 500-1 ■	500	0,1	В	9	
EMB 1200-1 =	1200	0,1	В		
EMB 2200-0 ▲	2200	1	В	•	
EMB 5.2K1 =	5200	1	В	9	
EMB 5.2K5 ▲	5200	5	В	•	

Overall dimensions WxDxH 170x240x39 mm. Net weight approx. 600 g.

KERN - Your DKD-Calibration Partner



Your KERN specialist dealer Distributed by:



ADVANCED APPLIED TECHNOLOGIES Contact Us:

Irl Ph: 01 4523432 UK Ph: 08452 30 40 30 Web: www.carlstuart.com Email: info@carlstuart.com