

LABORHITS

Bottle Carrying System DURAN® GL 45

Click and Carry

Safe Transport

The modular bottle carrying system fits standard laboratory bottles that have a GL 45 neck finish and is available in four sizes for 2, 5, 10 and 20 litres bottles. The carrying strap system is securely fastened around the bottle neck by the GL 45 bottle adapter. The adapter can only be removed by cutting.

Features

- Carrying handles allows easy and safe transport of bottles
- 5-litre bottle carrying system with anti-slip grip makes pouring liquids more comfortable and safer
- Fast and secure tightening of the polypropylene straps
- All materials are completely autoclavable
- Retrace Code for traceability and certification



DURAN® GL 45 Bottle Carrying System	colour	PK	Cat. No.
for 2 l bottles	blue	1	4.663 387
for 2 l bottles	purple	1	4.663 388
for 2 l bottles	yellow	1	4.663 389
for 2 l bottles	red	1	4.663 390
for 5 l bottles + anti-slip grip*	yellow	1	4.663 391
for 10 l bottles* + 2. handle	yellow	1	4.663 392
for 20 l bottles* + 2. handle	yellow	1	4.663 393

* incl. Silicone Bottle Base Protector

Bottles available separately.

Accessories Transport Dolly

Type	PK	Cat. No.
Metal dolly, for 10 l bottles, stainless steel	1	4.663 394
Metal dolly, for 20 l bottles, stainless steel	1	4.663 397
Handle for the 10 l metal dolly, stainless steel	1	4.663 395
Handle for the 20 l metal dolly, stainless steel	1	4.663 398
Silicone mat, yellow, for 10 l bottles	1	4.663 396
Silicone mat, yellow, for 20 l bottles	1	4.663 399

