


PDS No. 1882xx/ 188271-TRI	PRODUCT DATA SHEET				Page 1 of 1
Revision 12	Tube, 15 ml, PP, Screw Cap				 greiner bio-one
	Greiner Item-No. 1882xx /188271-TRI				
Valid for Item-No.:	188261 (sterile)	188262 (sterile)	188271 (sterile)	188272 (sterile)	188280 (sterile)
	188281 (sterile)	188283 (sterile)	188285 (sterile)	188271-TRI (sterile)	

1.	Description / Specification	
1.1	Description	PP tube, 15 ml, with screw cap, conical bottom, graduation and writing area, sterile 188271-TRI: triple packed
1.2	Dimensions	See Customer Drawing
1.3	Volume	Working volume: 14 ml
1.4	Material / Resin	<u>Tube</u> : PP (Polypropylene), free of heavy metal <u>Cap</u> : HDPE (High Density Polyethylene), free of heavy metal
1.5	Colour	<u>Tube</u> : 188261, -262, -271, -272, -271-TRI, -281, -285: natural ; print: graduation blue, writing area white 188280, -283: brown ; print: graduation blue, writing area white <u>Cap</u> : 188261, -262, -271, -272, -271-TRI, -280, -283: blue 188281, -285: white
1.6	Sterilization	188261, -271, -280, -281, -283, -285: SAL10 ⁻³ 188262, -272, -271-TRI: SAL10 ⁻⁶
1.7	Quality Control	- <u>Raw Material-Control</u> : physical testing - <u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	For single use only

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic
2.2	Temperature range	<u>Tube</u> : -196°C to +121°C <u>Cap</u> : -80°C to +100°C
2.3	Autoclavability	Yes (without cap)
2.4	Centrifugation, max. RCF	4000 x g: swinging-bucket rotor 15000 x g: fixed-angle rotor
2.5	Chemical Resistance	See homepage: https://www.gbo.com/en_INT/know-how-services/download-center.html
2.6	Shelf life	5 years after month of production
2.7	Other Information	-

3.	Packaging	188261, -262, -281, -283	188271, -272, -280	188285	188271-TRI
3.1	Pieces / bag or rack	50 in rack / bag	100 / bag	100 / bag	5 / bag (triple packed)
3.2	Pieces / box	500 pieces (10 racks)	1000	700	500
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)			
3.4	Other Information	Certificate of Quality			

4.	Other Information
4.1	188261, , -262, --271, -272, -271-TRI: Laboratory use only



Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 11	Date 21 August 2017	Date 22 August 2017	Date 22 August 2017	
Date 26.08.2015	Name S. Kaelberer	Name Dr. T. Schreiber	Name A. Schulz	

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.