

Optional for the 'HighTech' Series :



**Electronic module  
with analog connections**  
(order no. 8 900 100)

A with one input and two outputs for external programming or flow sensor or temperature recorder (current/voltage), scalable  
B Standby input  
C Alarm output

**HL + SL models (rear view)**



with connections for

- 1 External Pt1 00 sensor
- 2 RS232 / RS485 interface
- 3 Cooling machine / solenoid valve
- 4 **Electronic module with analog connections**
- 5 Connections for solenoid valve, etc.  
(230 V or 115 V)

**Online remote diagnosis**  
for all 'HighTech' circulators

Distributed by:



**Carl Stuart Limited**

ADVANCED APPLIED TECHNOLOGIES

**Contact Us:**

Irl Ph: 01 4523432

UK Ph: 08452 30 40 30

Web: [www.carlstuart.com](http://www.carlstuart.com)

Email: [info@carlstuart.com](mailto:info@carlstuart.com)