

# KU-3™ Viscometer

for Paints, Coatings and Inks

**ASTM D562 Compatible**  
(industry specification)

**NEW Magnetic Spindle Coupling**

**NEW User Interface**

**NEW Single Piece Pint and 1/2 Pint  
Can Adapter**

**Easy to use**  
no weights, simplifies an  
established test procedure

**LED Display Info:**

- Krebs Units
- Gram Units  
(Weight)
- Centipoise\*

Select Krebs or Grams or Centipoise

**Lock-In Test Results**  
with hold key

**Accuracy: ±1.0% of range**

**Repeatability: ±0.5%**

**Standard Krebs Spindle**

**Measurement range:**  
40 to 141 KU, 32 to 1099 gm, and 27  
to 5274 cP\*

**Adapter provided**  
for pint and half pint cans

**Accommodates quart cans**



## What's Included?

- Instrument
- Krebs-type Spindle (KU-1030)
- Adapter for Pint and Half-Pint Can

## Optional Accessories

- Paste Spindle (KU-1031)

## Applications

- Paints
- Coatings
- Adhesives
- Inks
- Pastes

\*Centipoise values based on the conversion from Krebs Units as defined in ASTM D562.

**Lab Unlimited**  
CARL STUART GROUP

Tallaght Business Park  
Whitestown, Dublin 24,  
Ireland  
D24 RFX3

Quatro House, Frimley Road,  
Camberley,  
United Kingdom  
GU16 7ER

Tel: (01) 4523432

Tel: 08452 30 40 30

Fax: (01) 4523967

Fax: 08452 30 50 30

E-mail: info@labunlimited.com

E-mail: info@labunlimited.co.uk

Web: www.labunlimited.com

Web: www.labunlimited.co.uk

ISO 9001:2008 CERTIFIED

